

Welcome to the PHC Webinar Series

This lecture on “**Next-Generation Sequencing for the Clinical Laboratory**” is given by
Karl V. Voelkerding, MD, FCAP



Your host is Jill Kaufman, PhD.
For comments about this webinar
or suggestions for upcoming
webinars, please contact
Jill Kaufman at jkaufma@cap.org

THE WEBINAR WILL BEGIN MOMENTARILY. ENJOY!

Karl Voelkerding, MD, FCAP



- Professor of Pathology at the University of Utah
- Medical Director for Genomics and Bioinformatics at the ARUP Laboratories
- Past President of the Association for Molecular Pathology
- Board certified in Clinical and Molecular Genetic Pathology
- His research interests include translation of nucleic acid based technologies into diagnostics with a current focus on complex genetic analyses by next generation sequencing



cap

Place sub-brand
here



Next-Generation Sequencing for the Clinical Laboratory

Karl V. Voelkerding, MD, FCAP

July 20, 2011

www.cap.org

v. #

Disclaimer

The College does not permit reproduction of any substantial portion of the material in this Webinar without its written authorization. The College hereby authorizes attendees of the CAP Webinar to use the pdf presentation solely for educational purposes within their own institutions. The College prohibits use of the material in the Webinar – and any unauthorized use of the College’s name or logo – in connection with promotional efforts by marketers of laboratory equipment, reagents, materials, or services.

Opinions expressed by the speaker are the speaker’s own and do not necessarily reflect an endorsement by CAP of any organizations, equipment, reagents, materials or services used by participating laboratories.

Disclosure

- I have nothing to disclose.

Outline

- **Progression: Gene Panels to Genomes**
- **Next Generation Sequencing Technology**
- **Bioinformatics**

Progression

**Whole
Genome**

**Whole
Exome**

**Multi-Gene
Diagnostics**

Increasing Complexity

The diagram illustrates a progression of genomic testing methods. At the bottom, a light gray wedge points to the right, labeled 'Increasing Complexity'. Above this wedge, four boxes are arranged in a staircase pattern from bottom-left to top-right. The first box, 'Multi-Gene Diagnostics', is light purple. The second, 'Whole Exome', is light green. The third, 'Whole Genome', is light gray. The final box at the top, 'Progression', is light green. All boxes have a dark blue border and a slight drop shadow.

Multi-Gene Diagnostics


Clinical Phenotype

Multiple Genes

Mutational Spectrum

Locus Heterogeneity

Allelic Heterogeneity



Multi-Gene Diagnostics

Cardiomyopathies

Hypertrophic

Dilated

Arrhythmias

10-35 Genes

Mitochondrial Disorders

Mitochondrial Genome

Nuclear Genes > 100 Genes

X-Linked Mental Retardation

~ 95 Genes

Multi-Gene Diagnostics

Hearing Loss

Retinopathies

Metabolic Disorders

Oncology

Structure/Function Complexes and Signaling Pathways

Multi-Gene Diagnostics

More Comprehensive Approach



Diagnosis



Prognosis



Treatment



Counseling

Multi-Gene Diagnostics

Technical Options

Sanger Sequencing of Individual Genes

Scanning and Sequencing

Multi-Gene Resequencing Microarrays

Next Generation Sequencing


Outline

- Progression: Gene Panels to Genomes
- **Next Generation Sequencing Technology**
- Bioinformatics

Sanger Sequencing


PCR followed by Cycle Sequencing with dNTPs/ddNTPs

Electrophoretic separation of chain termination products



Sequence DNA fragment library *in situ* in a flow cell

Massively parallel configuration




A T C G

[illegible]

Process

Genomic DNA or Enriched Target Genes



150-500
bp

End Repair and Adapter Ligation

+/-
PCR
"Fragment Library"

Adapter Fragment A Adapter

Adapter Fragment B Adapter

Adapter Fragment C Adapter


Process

“Fragment Library”




Clonal Amplification of Each Fragment

Emulsion Bead PCR




Surface Clusters



Sequencing of Clonal Amplicons in a Flow Cell


Process

Sequencing of Clonal Amplicons in a Flow Cell



Qualitative and Quantitative Information

Ref Seq



Multi-Gene Diagnostics

Genomic DNA

Enrichment

Target Genes

NGS Library Preparation

Next Generation Sequencing

Bioinformatics

Interpretation



Gene Enrichment Approaches

Genomic DNA

Amplification Based


Array Capture Based

PCR or LR-PCR
RainDance ePCR
Fluidigm
HaloGenomics

Solid Surface
or
In Solution

Enriched Genes

NGS



Gene Enrichment Approaches

Genomic DNA

Amplification Based

Array Capture Based

PCR or LR-PCR
RainDance ePCR
Fluidigm
HaloGenomics

Solid Surface
or
In Solution

Advantage: Enrichment Specificity

Advantage: Scalable to Exome

Drawbacks:
Not as Scalable
Instrument and Chip
Costs

Drawbacks:
Homologous Sequence
Capture
Manually Complex

Multi-Gene Diagnostics

Genomic DNA

Enrichment

Target Genes

NGS Library
Preparation



Next Generation Sequencing

Bioinformatics

Interpretation



NGS Library Preparation



Beckman SPRI-TE
1-10 Samples

Multi-Gene Diagnostics

Genomic DNA

Enrichment

Target Genes

NGS Library Preparation


Next Generation Sequencing

Bioinformatics


Interpretation



Next Generation Sequencing Technology



Illumina HiSeq 2000



Advantage: High Throughput

Drawbacks:
Batching
Sample Coordination



**Independent Flow Cells
8 Lanes per Flow Cell**

Multiple Panel Samples per Lane

1- 3 Exome(s) per Lane

1- 2 Genome per 8 Lanes

**New Platforms
Lower Throughput - Faster TAT
“Random Access”**



Illumina miSeq


Reversible Dye Terminators




Ion Torrent PGM

Monitors H⁺ Release


Phosphodiester Bond Formation



Multi-Gene Diagnostics




Multi-Gene Diagnostics – Parallel Testing



Re-Sequencing – Cardiomyopathy Genes

Variant g.34142190T>C in *TPM1*



Reference

LR-PCR
47%

Sanger

Progression

Whole Exome

**Multi-Gene
Diagnostics**

Increasing Complexity



The diagram illustrates a progression of genomic testing methods. It features three rectangular boxes arranged in a staircase pattern from bottom-left to top-right. The bottom-left box is purple and labeled 'Multi-Gene Diagnostics'. The middle box is green and labeled 'Whole Exome'. The top-right box is light green and labeled 'Progression'. Below these boxes is a large, light gray wedge that tapers from left to right, with the text 'Increasing Complexity' centered within it. All elements have a dark blue border.

Human Exome




“Journey to the Center of the Genome”

~ 30 Megabases (~ 1% of the genome)

~ 180,000 exons (~ 20,500 genes)

Harbors “Majority” of Mendelian Mutations

Gene Discovery
~ 40 Publications
July 2011



Genomic DNA



Library Preparation

Next Generation Sequencing Library



Hybridize to Exome Capture Probes


Exome Enriched Library



**Next Generation
Sequencing**

Bioinformatics Analysis

Exome Sequencing – Cardiomyopathy Genes



Variant g.34142190T>C in *TPM1*

Reference

LR-PCR
47%


Exome
47%

Sanger

GAGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
GCCAAGCACATTGCTGAAGATGCCGACCGCAAGACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAATGCGGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
GCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC
CCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACC

GAGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
GAGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGG
GAGGCCAAGCACATCGCTGAGGATTACGACCGCAAATATGAAGAGGT
GAGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGAC
AGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACA
CCAAGCACATTGCTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCC
CAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCA
CTGAAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
TGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
AAGATGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
TGCCGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
CGACCGCAAATACGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG
GACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG

GAGGCCAAGCACATTGCTGAAGATGCCGACCGCAAATATGAAGAGGTCAGATCCTGGGGCCCCAAAGCCTTGTGGACACCCAG



Progression

Whole Genome

Whole Exome

**Multi-Gene
Diagnostics**

Increasing Complexity

The diagram illustrates the progression of genomic testing. It features four rectangular boxes arranged in a staircase pattern from bottom-left to top-right. The first box is light blue and labeled 'Multi-Gene Diagnostics'. The second box is light green and labeled 'Whole Exome'. The third box is light green and labeled 'Whole Genome'. The fourth box is light green and labeled 'Progression'. At the bottom, a light gray wedge with a blue outline points from left to right, labeled 'Increasing Complexity'.

Genomic DNA




Fragmentation



QC




Covaris
Acoustic Wave



Library Preparation - Illumina

Fragment End Repair

Adapter Ligation




Beckman SPRI-TE

PCR



QC

Process



BioAnalyzer

PCR Amplified Library




Gel Electrophoresis




Size Selection



Gel Purification



QC



Gel Purified Library



Quantitative PCR




Dilution/Denaturation




5-7pM

Flow Cell Cluster Generation



Illumina cBot

Process



cBot - Top View

Sequencing HiSeq 2000

Process

↓
“First Base” Report

QC

Cluster Densities/Intensities


↓
Real Time Analysis - Cycle 25+

QC

2 X 100bp
Pair End


→ Target
~ 100+ Gb

△
→ FastQ Files



Whole Genome Sequencing


Chr 10: g.43,615,633C>G in *RET*



Whole Genome Sequencing



Chr X: 3bp deletion in *FOXP3*

Mother



p-19


Son




Outline

- Progression: Gene Panels to Genomes
- Next Generation Sequencing Technology
- **Bioinformatics**

Next Generation Sequencing Bioinformatics



Alignment Considerations



Alignment Softwares


Academic Softwares

- Command Line
- Free
- Genome community support

```
maq assemble [-sp] [-m maxmis] [-Q maxerr] [-r  
hetrate] [-t coef] [-q minQ] [-N nHap]  
out.cns in.ref.bfa in.aln.map 2> out.cns.log
```

Commercial Softwares

- Feature rich user interface
- \$\$\$
- Company support



Alignment Softwares

Commercial Softwares




Academic Softwares

[BFAST](#) (UCLA)

[BWA](#) (Sanger Institute)

[SAMtools](#) (Sanger Institute)

Next Generation Sequencing Bioinformatics




Variant Annotation

- ✓ dbSNP ← 1,000 Genome Project
- ✓ OMIM
- ✓ Human Genome Mutation Database
- ✓ Locus Specific Databases
- ✓ Literature and Internet
- ✓ Functional Prediction Programs

PolyPhen

SIFT

Convergence of Chemistry + Bioinformatics



Exome Sequencing

Variant Calling
(15-20,000)

Filter Out
Common
Variants
(750-1000)

Genes/Regions
Family/SNP
Arrays (10-200)

Variant
Annotation

Candidate
Genes

Approach 1

Nonsense
Splicing &
Frame shifts

Missense:
Protein Function
Predictions


Functional
Relevance


Approach 2

Evolutionary
Constraint

Approach 3

VAAST
(Variant Annotation,
Analysis, and
Selection Tool)





Courtesy
Mark Yandell
Martin Reese

VAAST – Probabilistic Candidate Gene Finder

**Allele Frequencies
Cases and Controls**

**AA Substitution Analysis
Model Variant Severity**




**Combined Likelihood
Framework**




**Identify Aberrant Variant
Combinations
Compromise Gene Function**

VAAST – Probabilistic Candidate Gene Finder



VAAST – Miller Kindred Quartet



Miller + PCD



= Affected

Exome Sequence Data

Two Siblings



VAAST

VAAST – Miller Kindred Quartet

Caucasian Only (65 genomes)				
Genome-wide Significant Genes	DHODH		DNAH5	
	Rank	P-Value	Rank	P-Value
17	14	9.93E-07	19	5.79E-05


Mixed Ethnicities (189 Genomes)				
Genome-wide Significant Genes	DHODH		DNAH5	
	Rank	P-Value	Rank	P-Value
9	4	7.60E-09	5	1.18E-08

Miller : *DHOD*

PCD: *DNAH5*

GWS Alpha is 2.4
E-6

VAAST – Miller Kindred Quartet



Miller + PCD



= Affected

Exome Sequence Data

Two Siblings

Parents

VAAST


VAAST – Miller Kindred Quartet

One



Two

Genomics Clinical Research Program



Summary

- **Progression: Gene Panels to Genomes**
- **Next Generation Sequencing Technology**
- **Bioinformatics**

Acknowledgements

Genomics-Bioinformatics

**Rebecca Margraf Jacob Durtschi
Emily Coonrod Perry Ridge**

U of Utah Genetics

Mark Yandell Lynn Jorde

voelkek@aruplab.com

Omicia

Martin Reese

**ARUP Laboratories Institute
for
Clinical and Experimental Pathology**

Next in the Series of Free PHC Webinars

- **How to Have Successful Patient Interactions, Wednesday, August 17th, 11:00-12:00 pm CT**
 - Mary Ann Abrams, MD, MPH & Barbara Savage, MT(ASCP)
- Go to www.cap.org/institute For All Upcoming Webinars!
- Past Webinars Available Now Online at www.cap.org/institute
 - Accountable Care Organizations
 - Whole Genome Analysis as a Universal Diagnostic
 - How to Build and Fund a Financially Viable Molecular Lab
 - Cancer: The Critical Role of Pathology
 - Molecular Markers in Breast Cancer
 - Bethesda System: Integrating Cytology and HPV Molecular Testing
 - Molecular Diagnosis for Lung Cancer Patients
 - Molecular Diagnosis for Colorectal Cancer Patients

CAP Events of Interest

- Don't Forget to Register for **CAP'11 – THE Pathologists' Meeting** – September 11 – 14, 2011 held at the Gaylord Texan in Grapevine, Texas!
 - Go to www.cap.org/CAP11 or call 1-800-967-4548. International attendees please call 1-847-996-5891.

For more information go to www.cap.org/CAP11

Tuesday, Sept 12th:

TP120 Breakfast Workshop – Hot Topics in Pathology: What Every Community Pathologist Should Know About Clinical Requests for Molecular Tests (6:30-7:45 am)

Faculty--Samuel K. Caughron, MD, FCAP

Frederick L. Kiechle, MD, PhD, FCAP

Michael S. Brown, MD, FCAP

ST109 Companion Diagnostics for Targeted Therapy in Cancer (2:00-5:30 pm)

Faculty--Sanja Dacic, MD, PhD, FCAP

David Hicks, MD, FCAP

Jeffrey Kant, MD, PhD, FCAP

Wednesday, Sept 13th:

ST110 Direct-to-Consumer Genetic Testing: Staying Ahead of Patients in This Current Trend (8:00-9:00 am)

Faculty--Nazneen Aziz, PhD

Elizabeth A. Mansfield, PhD

ST111 What's in It for Me? Using Technology to Become a Diagnostic Hero (8:00-11:30 am)

Faculty--Kenneth J. Bloom, MD, FCAP

John W. Turner, MD, FCAP